

Transducer Disinfection



Terason's transducers

Our family of lightweight and ergonomic transducers is ideal for easy scanning


15L4 Smart Mark™
Linear Array



16L5 Smart Mark™
Linear Array



8L2 Smart Mark™
Linear Array



4V2 Phased Array



8V3 Phased Array



5C2 Curved Array




9MC3 Curved Array



8EC4 Endocavity



16HL7 High Frequency
Linear Array



Pedoff




8TE3 Transesophageal



10EC4 Endocavity



8BP4 Endocavity



Transducer Disinfection

The following products are approved for use with Terason uSmart transducers:

Ultrasound Gels	Disinfectants	High-Level Disinfectants*
<ul style="list-style-type: none"> Aquagel Aquasonic Clear Clear Image Clear wave Scan Ultrapasonic 	<ul style="list-style-type: none"> Cavi Wipes Transeptic Spray Sono Wipes P.I. Spray II Protex Wipes Sani-Coth Wipes Sani-Cloth HB Wipes Sani-Cloth Plus Wipes 	<ul style="list-style-type: none"> Cidex Cidex OPA Cidex Plus Metricide 14 Metricide 28 Wavicide <p>*For use only on endocavity and TEE transducers</p>

Trophon EPR system has been approved for the following:

15L4A, 16L5, 5C2A, 8V3, 9MC3, 8EC4A, 15L4, 4V2A, 8V3A, 8L2

Trophon2 system has been approved for all probes listed for the Trophon EPR system, and adds 16HL7



Terason Division • Teratech Corporation
77 Terrace Hall Avenue, Burlington, MA 01803 • 1-866-TERASON • 781-270-4143 • terason.com

terason.com